


SHL-US20

NEBULIZATOR ULTRASONIC control total al respiratiei

- Distribuire rapida si eficienta a medicatiei
- Compact si ultraportabil
- KIT complet de nebulizare
- 2 Moduri de functionare: - Mod NORMAL
- Mod INTENS
- Reglarea volumului de aer (control complet al tratamentului)
- Capacitate mare a cupei de medicatie
- Oprire automata
- Mod de alimentare: - alimentator AC la priza (inclus)
- Usor de intrebuintat
- LED-uri de semnalizare



NEBULIZATOR ULTRASONIC SHL-US20



CARACTERISTICI

- Distribuire rapida si eficienta a medicatiei
- Compact si ultraportabil
- KIT complet de nebulizare COMPUS din:
 - cupe medicatie 5 buc.
 - pipeta
 - piesa de gura
 - masca pentru adulti
 - masca pentru copil
 - adaptor pentru masca
 - alimentator pentru alimentare la priza
 - husa pentru transport
- 2 Moduri de functionare: - Mod NORMAL (20 min.)
 - Mod INTENS
- Reglarea volumului de aer pentru control complet al tratamentului

SPECIFICATII TEHNICE

- Diametrul mediu al particulelor nebulizate (MMAD): 1-5 μm
- Capacitate maxima a cupei de medicatie: max. 8 ml
- Rata medie de nebulizare: 0.25ml/min – 0.7 ml/min (functie de densitatea medicamentului)
- Oprire automata: dupa aprox. 3 min. de neutilizare
- Dimensiuni: 63 mm x 45 mm x 152 mm
- Greutate: aprox. 120 gr.
- Frecventa ultrasunetelor: 2.5 MHz
- Nivel de zgomot: ≤ 55 dBA
- Lichid rezidual: max. 1 ml
- Tensiune de alimentare: 100-240V, 50/60 Hz, 500 mA

SOLUTII PENTRU O SANANTATE PERFECTA

